nf-core/ x scrnaseq


Chris Hakkaart - Seqera

nf-core/scrnaseq


nf-core/scrnaseq is a bioinformatics best-practice analysis pipeline for processing 10x Genomics single-cell RNA-seq data.

This is a community effort in building a pipeline capable to support:

- Alevin-Fry + AlevinQC
- STARSolo
- Kallisto + BUStools
- Cellranger
- UniverSC







--outdir results

-profile test,docker \

nextflow run nf-core/scrnaseq \

```
params {
   config profile name = 'Test profile'
   config profile description = 'Minimal test dataset to check pipeline function'
   max_cpus = 2
   max memory = '6.GB'
   max time = '6.h'
   input = 'samplesheet-2-0.csv'
   fasta = 'GRCm38.p6.genome.chr19.fa'
   gtf = 'gencode.vM19.annotation.chr19.gtf'
   aligner = 'star'
   protocol = '10XV2'
   validationSchemaIgnoreParams = 'genomes'
process {
   withName: '.*:CELLRANGER COUNT' {
```

maxForks = 1

```
params {
   config_profile_name = 'Full Test profile'
   config profile_description = 'Full test dataset to check pipeline function'
   input
               = 'samplesheet_2.0_full.csv'
   // Genome references
               = 'GRCh38'
   genome
   aligner = 'star'
   protocol = '10XV2'
   validationSchemaIgnoreParams = 'genomes'
```

"5000"

samplesheet.csv

sample	fastq_1	fastq_2	expected_cells	seq_center
Sample_X	Sample_X_S1_L1_R1.fastq.gz	Sample_X_S1_L1_R2.fastq.gz	"5000 "	"Site 1"
Sample_Y	Sample_Y_S1_L1_R1.fastq.gz	Sample_Y_S1_L1_R2.fastq.gz	"5000"	"Site 2"
Sample_Y	Sample_Y_S1_L2_R1.fastq.gz	Sample_Y_S1_L2_R2.fastq.gz	"5000"	"Site 3"

```
$ nextflow run scrnaseq/ -profile test,docker --outdir results
NEXTFLOW ~ version 23.10.0
Launching `scrnaseq/main.nf` [ecstatic miescher] DSL2 - revision: bd078e3da5
executor > local (14)
[e9/3dc2c1] process > NFCORE SCRNASEO:SCRNASEO:INPUT CHECK:SAMPLESHEET CHECK (samplesheet-2-0.csv)
                                                                                                       [100%] 1 of 1 🗸
[42/930a0d] process > NFCORE SCRNASEQ:SCRNASEQ:FASTQC CHECK:FASTQC (Sample Y)
                                                                                                       [100%] 2 of 2 🗸
                                                                                                       [100%] 1 of 1 🗸
[cc/df9dca] process > NFCORE SCRNASEQ:SCRNASEQ:GTF GENE FILTER (GRCm38.p6.genome.chr19.fa)
[81/d73a70] process > NFCORE SCRNASEQ:SCRNASEQ:STARSOLO:STAR GENOMEGENERATE (GRCm38.p6.genome.chr19.fa) [100%] 1 of 1 🗸
[8c/2c18e3] process > NFCORE SCRNASEO:SCRNASEO:STARSOLO:STAR ALIGN (Sample Y)
                                                                                                       [100%] 2 of 2 🗸
[0b/1ca671] process > NFCORE SCRNASEQ:SCRNASEQ:MTX CONVERSION:MTX TO H5AD (Sample X)
                                                                                                       [100%] 2 of 2 🗸
[cb/8a97ee] process > NFCORE SCRNASEO:SCRNASEO:MTX CONVERSION:CONCAT H5AD
                                                                                                       [100%] 1 of 1 🗸
[ae/e8c1c4] process > NFCORE SCRNASEQ:SCRNASEQ:MTX CONVERSION:MTX TO SEURAT (Sample Y)
                                                                                                       [100%] 2 of 2 🗸
[b9/35129b] process > NFCORE SCRNASEO:SCRNASEO:CUSTOM DUMPSOFTWAREVERSIONS (1)
                                                                                                       [100%] 1 of 1 🗸
[ef/2862c4] process > NFCORE SCRNASEO:SCRNASEO:MULTIOC
                                                                                                       [100%] 1 of 1 🗸
-[nf-core/scrnaseq] Pipeline completed successfully-
Duration : 7m 17s
CPU hours : 0.3
```